

## NEWS RELEASE 6/15/2000

### NEW PRODUCT RELEASE

## **ELECTROPHORESIS POWER SUPPLIES ARE PROGRAMMABLE AND COME WITH A 3 YEAR WARRANTY**

Topac introduces the Consort range of power supplies for electrophoresis. This range of power supplies covers all the common applications in electrophoresis including DNA sequencing, vertical and horizontal gels, blotting, isoelectric focusing and electro elution. Each power supply can function in constant current, constant voltage or constant power mode and can be programmed with up to 9 steps. The uncluttered front panel, membrane keypad, and clear digital display make the power supplies easy to configure and use. Each unit provides 4 sets of output terminal so that four individual apparatus (sequencers, blotters etc.) can be run simultaneously.



Voltage ranges from 150 Volts to 6000 volts are available. The power supplies also feature very durable construction, a high level of safety engineering and a three-year warranty.

For more information visit the web page for this product at

[http://www.topac.com/ep\\_es\\_E800.html](http://www.topac.com/ep_es_E800.html)

For a digital copy of this text and images visit

<http://www.topac.com/pressinfo.html>

Or contact Antoni Drybanski 781 740 8778 email [toni@topac.com](mailto:toni@topac.com)